



Vial Adaptor

With Toxi-Guard® technology

TECHNICAL DATASHEET

The Chemfort™ Vial Adaptor family is PVC-free and offers 4 sizes that firmly connect onto any standard vials, that have a cap dimension of 13 & 20 mm (ref. MG245248), 28 mm (ref. MG245266) or 32 mm (ref. MG245269).

Equipped with the Toxi-Guard® technology, an immediate pressure equilibration is ensured, reducing the extra step of pre-priming.

The Toxi-Guard® system contains a 100% activated carbon drug binding matrix (charcoal) and a 0.2 micron hydrophobic and oleophobic membrane. The membrane prevents airborne contaminants and particles from entering the drug or fluid path. The active carbon prevents the escape of cytotoxic and other hazardous drug vapors, aerosols and droplets into the environment during drug reconstitution.

Once connected to the Chemfort™ Syringe Adaptor/Syringe Adaptor Lock, a safe drug transfer is made possible from the vial to another container (e.g. a semi-rigid bottle, an infusion bag, an elastomeric pump or a medication cassette reservoir for an ambulatory infusion pump).

The intuitive design makes it easy to use: with minimal effort both Chemfort™ components can be connected and a "click" sound will indicate a safe connection. A pinch on the SA/SAL wings will disconnect the two components. The self-sealing elastomers prevent drug leakage and droplet formation on the outer septa surfaces. Also, in case of multiple-use vials, a protective cap can be used for safe storage of the vial.

The Chemfort™ portfolio has no known incompatibilities with hazardous drug solvents.

Indications for Use:

Chemfort™ is a Closed System Transfer Device (CSTD) that mechanically prohibits the release of drugs, including antineoplastic and hazardous drugs, in vapor, aerosol, or liquid form during preparation, reconstitution, compounding and administration, minimizing exposure of individuals, healthcare personnel and the environment to hazardous drugs.

Chemfort™ prevents the introduction of microbial and airborne contaminants into the drug or fluid path for up to 7 days.

TECHNICAL DATASHEET

Order Information

| | Vial Adaptor 20 mm with 13 mm Converter | Vial Adaptor 28 mm | Vial Adaptor 32 mm |
|----------------------------|---|-----------------------|-----------------------|
| Article number / Reference | MG245248 | MG245266 | MG245269 |
| User box quantity | 50 pcs | 20 pcs | 25 pcs |
| Transport box quantity | 200 pcs | 80 pcs | 100 pcs |

Product Characteristics

| Product technical data | | | |
|-------------------------------|---|---|---|
| Cap dimension of fitting vial | 13 & 20 mm | 28 mm | 32 mm |
| External fitting | Tight circular grip around the neck of the vial | Tight circular grip around the neck of the vial | Tight circular grip around the neck of the vial |
| Spike type | 2-way spike with separate air and fluid path | 2-way spike with separate air and fluid path | 2-way spike with separate air and fluid path |
| Protective cap | Yes (white) | Yes (white) | Yes (blue) |
| Fluid path volume | 0.15 ml | 0.15 ml | 0.15 ml |
| Membrane pore size | 0.2 micron | 0.2 micron | 0.2 micron |
| Septum piercings | Up to ten times | Up to ten times | Up to ten times |
| Shelf life | 3 years | 3 years | 3 years |
| Sterilization | Ethylene Oxide | Ethylene Oxide | Ethylene Oxide |

Tolerates disinfection with alcohol 70% and chlorhexidine in alcohol (5 mg/70% alcohol)

Material Specifications

| Component | Material | Material | Material |
|-----------------------|--------------------|--------------------|--------------------|
| Spike | Polypropylene | Polypropylene | Polypropylene |
| Injection point cover | Polyethylene | Polyethylene | Polyethylene |
| Elastomeric sealing | Polyisoprene | Polyisoprene | Polyisoprene |
| 13 mm convertor | PET | | |
| Protective cap | Polypropylene | Polypropylene | Polypropylene |
| 0.2 micron membrane | Acrylic co-polymer | Acrylic co-polymer | Acrylic co-polymer |
| Charcoal filter | Charcoal Zorflex® | Charcoal Zorflex® | Charcoal Zorflex® |
| Filter base | Polypropylene | Polypropylene | Polypropylene |
| Filter cover | PET | PET | PET |

Free of PVC, DEHP and Latex. Non-Pyrogenic. Free from lead, cadmium, arsenic and chrome. Meets standards and regulations with regards to the environment

Regulation

| |
|--|
| CE 0483 |
| FDA 510(k) Cleared under ONB product code |
| Disposal of components should be according to local regulations and procedures for handling of hazardous drugs |

Symbols taken from ISO 15223-1:2016



Sterilized using ethylene oxide



Non-pyrogenic



Do not re-sterilize



Do not re-use



Do not use if package is damaged



Keep dry



Keep away from sunlight



Donawa Lifescience Consulting Srl
Piazza Albania 10,
00153 Rome, Italy

CE 0483

SIMPLIVIA™

Simplivia Healthcare Ltd.
North Industrial Zone
P.O. Box 888
Kiryat Shmona, Israel 1101801



www.simplivia.com



info@simplivia.com



Follow us



SIMPLIVIA